

## 2 nozzles

# processing 0

HX3042/00

# A powerful stream of water

For a regular clean

Use the standard nozzle for a regular clean. It uses a single stream of water to remove debris and clean thoroughly between your teeth.

## An effortless clean with advanced technology

- F1 Standard nozzle for a regular clean
- Guidance tip ensures correct placement

**PHILIPS** 



## **Specifications**

#### Items included

• F1 Standard nozzle: 2

#### Ease of use

- · Nozzle attachment: Easily snaps on and off
- Nozzle replacement: Every 6 months, for good

#### hygiene

#### **Compatibility**

· Compatible with: any Philips Sonicare Power

## Highlights

A powerful stream of water
Use the standard nozzle for a regular clean. It uses a single stream of water to remove debris and clean thoroughly between your teeth.

### **Guidance tip**

The slim angled nozzle and guidance tip make finding the correct placement easy. Position the nozzle tip just above the gum line, press gently so that the nozzle tip is in contact with the gum line and the teeth, and slide from one tooth to the next.



Issue date 2021-07-11

© 2021 Koninklijke Philips N.V.

All Rights reserved.

Version: 4.0.1

Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V.

or their respective owners.

EAN: 08 71010 39470 28 www.philips.com